

IMMUNOLOGY: MANCINI TECHNIQUE

Reference : MANCINISUB



Substitute Mancini immunology kit for 20 pairs

L'immunité adaptative

Illustration of MANCINI technique with substitution products.
The results are obtained during the PW session.

Composition of the kit (12 to 20 pairs of students):
20 test dishes and 10 training dishes

- 30 petri dishes - 5.5 cm diameter
 - 30 die cutter droppers
 - 2 g of agar
 - 2 tubes containing the substitute products
 - Technical and educational instructions
- Necessary material:
- Variable volume micropipette 1-10 μL
 - 10 μL micropipette cones

- 250 ml Erlenmeyer
- Microwave or hotplate
- Humid chamber
- Graduated ruler for measurement