

INTRODUCTION TO MICROBIOLOGY

Reference : C/INIMIC



Introduce students to the cultivation of a bacterial strain.
We selected the strain *Halobacterium salinarum* because it has the advantage of growing on a very salty culture medium which prevents the growth of any other microorganism.
Students therefore handle it safely.

Composition :

- 1 strain of *Halobacterium salinarum* on a Petri dish
- 20 dishes of medium, diameter 90mm
- 20 sterile inoculating loops
Technical and educational instructions available on our website.

Conservation: 1 month at 4° C

Number of pairs of students: for 20

Necessary material: 30° C Oven , Bunsen burner (sterile handling)