


KIT PROTÉINES POUR LA CUVE VERTICALE

PRODUCT SHEET

Reference : KITPROT



Protein kit for the vertical tank

Protein molecular weights are determined by vertical polyacrylamide gel analysis. The molecular weight of pre-colored proteins, initially unknown, can be determined by comparison with the size marker corresponding to standard colored proteins.

Composition :

Instructions, denatured proteins, standard protein markers, gel loading solution, buffer.

Electrophoresis 60 min. Staining 20 min*.

Optional* decolorization 15 to 60 minutes.

*reading possible without staining