

MELTING POINT DEVICE

Reference : SMP11



Easy to use

With a manually adjustable heating rate, the SMP11 quickly heats up samples to the melting temperature. Up to 3 samples can be examined simultaneously.

Supplied with a mercury-free thermometer, there are no health risks in case of breakage.

The capillary tubes containing the samples to be tested are placed in the apparatus. A magnifying glass and a lighting system allow the visualisation of the sample and the reading of the temperature.

Delivered with 1 thermometer without mercury and a batch of 100 capillary tubes.

Temperature range: ambient to 250 ° C

Temperature rise: 1 to 10 ° C / min

Dimensions / Weight: 370 x 140 x 110 mm / 1.7 kg