


OBSERVATION OF MITOSES OF ONION

Reference : C/MITOIGN



Have students perform a microscopic preparation and observe mitotic figures from onion epidermis.

Composition:

- 1 dropper flask of 10ml of a sample preparation solution
- 1 dropper flask of 10ml of sample staining solution
- 1 dropper flask of 10ml of a solution for observation of the sample
Educational and technical instructions available on our website.

Conservation: 6 months at room temperature

Number of pairs of students: 20

Necessary material:

- Onion, garlic or shallot
- Beaker
- Watch glass
- Blade, cover slip

- Magnification microscope x 10, x 40
- Scalpels
- Scissors
- Fine forceps