

References : IMVD, IMVS



## The true Ouchterlony L'immunité adaptative

"Formation of an antibody-antigen complex"

Cognitive objectives:

Perform an Ouchterlony test.

To highlight the specificity of the humoral mediated immune reaction by showing the formation of antigen-antibody complex.

Benefits:

The preparation of the kit is easy and fast, just boil the medium in the microwave and then distribute it in Petri dishes.

The result is guaranteed, it is observed after 24 hours and can be kept for 2 weeks in the refrigerator.

Conservation:

6 months in the refrigerator

Note: The serums keep very well if frozen between two actions.

Double version (for 40 pairs of students)

**Réf : IMVD**

Simple version (for 20 pairs of students)

**Réf : IMVS**