

REVEALS GAMMA-GLOBULINS

Reference : 840582



This kit can detect the presence of blood on fabric or other materials by performing a Kastle-Meyer test. A wet cotton swab is rubbed on the tissue to be analysed and the reagents are poured over it. The appearance of colour indicates the presence of human or animal blood. Positive control strips are used to verify that the technique is working and to obtain the control colour.

Composition of the kit::

- 50ml of phenolphthalein
 - 50ml of ethanol (90%)
 - 50ml of hydrogen peroxide (3%)
 - 5 positive control strips simulating blood-soaked tissue *
 - 50 cotton swabs
 - Technical and educational instructions
- * They let you take a sample with the cotton swab and give a positive test