

SUBSTITUTE OUCHTERLONY

Reference : IMSUB



L'immunité adaptative

This kit makes it possible to simulate an Ouchterlony test and to obtain a result that is clear and identical to a real test in less than an hour.

Composition for 30 pairs of students:

- 30 Petri dishes
 - 10 glass die cutters
 - 10 calibrated droppers
 - 1 empty bottle
 - 1 bag of agar
 - 2 tubes of solution mimicking rabbit serum immunised against BSA (NaOH)
 - 2 tubes of solution mimicking the BSA (ZnSO₄)
- Technical and educational instructions available on our website.

Necessary material:

- Box and paper towel to form a humid chamber

- Distilled water
- Microwave or water bath and anti-heat glove
- Microtubes
Conservation: 1 year
Room temperature