

Reference : PS-2179



This colorimeter is specially designed for water analysis. It makes it possible to very easily and accurately measure the ionic concentration of aqueous samples for a choice of 8 factors. It works in combination with the reagent kits opposite.

It is very simple to use:

The student has reagent ampoules (corresponding to the ion to be titrated), just break the tip of the ampoule with a safe system (provided) and then place the ampoule in the sample to be analysed.

The bulb fills automatically, then all you have to do is shake it to mix and place the bulb in the colorimeter. The measurements are displayed in seconds.