CRYO BEADS

Reference : CRYOBILLES





Cryo beads are a simple way to keep important microorganisms free of contamination. The product is composed of 25 porous beads in a vial filled with a broth containing glycerol. Microorganisms stick to the porous surface of the beads. Excess broth is aspirated, the flasks can then be stored at -70 $^{\circ}$ C in a freezer.



15 Av. des grenots - 91150 Etampes Tel. : +33 (0)1 69 92 26 72 - Fax : +33 (0)1 69 92 26 74 - Email : info@sordalab.com