

INSTANT DEHYDRATED MEDIUM

References : MDRBE4L, MDRBE1L, MDR1L, MDR4L



No preparation

The preparation of this medium requires no cooking, just add cold water. A tube of culture medium can thus be prepared in less than one minute.

Consisting of all the components necessary for the proper development of *Drosophila* and an antifungal, very resistant to dehydration, this environment will facilitate the establishment of all your cultures and crosses. 1 dose of one litre allows the preparation of 70-75 culture tubes.

This medium is available in natural colour or with a blue dye that makes it easier to detect the larvae and arouses students' curiosity.