
MEDIA



References: MEUG, MEUG.00



Bottles of sterile liquid medium ready for use

Three types of medium are offered:

For the culture of euglena (MEUG) and chlorella (MCLO)The MEUG culture medium contains lactic acid, it is used to obtain a rapid growth of euglena. (The euglena breathe at the beginning of growth then photosynthesise at the end of the culture)

The MEUG.00 culture medium does not contain lactic acid (euglena photosynthesise).

Euglena medium + lactic acid , sterile

Réf : MEUG

Euglena medium without lactic acid, sterile

Réf: MEUG.00

