

# RACK OF CONES

References : RACKB, RACKJ, RACKG



Rack of 10-200  $\mu$ l yellow cones  
Cones presented in a 96-point rack.  
Polypropylene rack with transparent polycarbonate lid.  
Autoclavable at 110° C.

**Réf : RACKJ**

Rack of 100-1 000  $\mu$ l blue cones  
Cones presented in a 96-point rack.  
Polypropylene rack with transparent polycarbonate lid.  
Autoclavable at 110° C.

**Réf : RACKB**

Rack of 2-20  $\mu$ l gray cones  
2-20  $\mu$ L Cones presented in a 96-point rack.  
Polypropylene rack with transparent polycarbonate lid.  
Autoclavable at 110 ° C. The rack can also be placed in the Tip SystemBox.  
You can fill it with cones up to 20 $\mu$ L after using up those delivered.

**Réf : RACKG**